

Applications Note

GPR series AUTOMATIC REFRACTOMETER

The following substances have been measured successfully on the GPR Automatic Refractometer:

ADHESIVES - PVA & hot melt

ESSENTIAL OILS, FLAVOURS & FRAGRANCES

PLASTICISERS

OILS

PHARMACEUTICALS e.g. Penicillin

WAXES & POLISHES

FOODSTUFFS - Gelatin
 Milk & milk products *
 Sucrose/Glucose
 Tomato juice, paste & puree
 Mincemeat *
 Rice Pudding *
 Lemon Curd *
 Jam *
 Tartar sauce *
 Brown sauce
 Peanut butter *
 Honey (set or cloudy type) *
 Malt extract
 Caramel *
 Fruit juice

** on Emulsion Scale setting*